# **SALTER**

## **PERSONAL SCALES**

### **Instructions and Guarantee**

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**PERSONENWAAGE**

**BALANZA PERSONAL**

**BILANCIA PESAPERSONA**

**BILANÇÀ PESAPERSONA**

**BADEROMSVEKT**

**PERSONENWEEGSCHAAL**

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**PERSONVÅG**

**PERSONVÆGT**

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**OSOBNÍ VÁHA**

**KİŞİSEL TERAZİ**

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## INSTRUCTIONS FOR USE

Your Salter Scale is a quality item which for best results should be used in the following way:

1. Use on a firm flat surface.
2. Set scale to zero using the serrated wheel.
3. Always weigh at the same time each day, without clothes or any footwear.
4. Stand as still as possible with feet evenly spaced on the scale.
5. To clean, only use a damp cloth. NEVER use chemicals or abrasives.
6. For maximum durability, keep scale dry and avoid direct heat or sunlight.

## GUARANTEE

This product is intended for domestic use only. Salter will repair or replace the product, or any part of this product, (excluding batteries) free of charge if within 15 years of the date of purchase, it can be shown to have failed through defective workmanship or materials. This guarantee covers working parts that affect the function of the scale. It does not cover cosmetic deterioration caused by fair wear and tear or damage caused by accident or misuse. Opening or taking apart the scale or its components will void the guarantee. Claims under guarantee must be supported by proof of purchase and be returned carriage paid to Salter (or local Salter appointed agent if outside the UK). Care should be taken in packing the scale so that it is not damaged while in transit. This undertaking is in addition to a consumer's statutory rights and does not affect those rights in any way. For UK Sales and Service contact HoMedics Group Ltd, PO Box 460, Tonbridge, Kent, TN9 9EW, UK. Helpline Tel No: (01732) 360783. Outside the UK contact your local Salter appointed agent.

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