SALTER

Body Analyser & Scale

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HOW DOES THIS SALTER SCALE WORK?

This Salter Scale uses BIA (Bio Impedance Analysis) technology which passes a tiny electrical impulse through the body to determine fat from lean tissue, the electrical impulse cannot be felt and is perfectly safe. Contact with the body is made via stainless steel pads on the platform of the scale.

This method simultaneously calculates your personal weight, body fat, total body water, muscle mass and BMI giving you a more accurate reading of your overall health and fitness.

This scale stores the personal data of up to 10 users. As well as being an analyser scale, this scale can be used as a conventional scale.

BODY FAT - WHAT DOES IT MEAN?

The human body is made up of, amongst other things, a percentage of fat. This is vital for a healthy, functioning body - it cushions joints and protects vital organs, helps regulate body temperature, stores vitamins and helps the body sustain itself when food is scarce. However, too much body fat or indeed too little body fat can be damaging to your health. It is difficult to gauge how much body fat we have in our bodies simply by looking at ourselves in the mirror.

This is why it is important to measure and monitor your body fat percentage. Body fat percentage gives you a better measure of fitness than weight alone – the composition of your weight loss could mean you are losing muscle mass rather than fat - you could still have a high percentage of fat even when a scale indicates 'normal weight'.

BODY WATER - WHY MEASURE IT?

Body water is the single most important component of body weight. It represents over half of your total weight and almost two thirds of your lean body mass (predominantly muscle).

Water performs a number of important roles in the body:

All the cells in the body, whether in the skin, glands, muscles, brain or anywhere else, can only function properly if they have enough water. Water also plays a vital part in regulating the body's temperature balance, particularly through perspiration.

The combination of your weight and fat measurement could appear to be 'normal' but your body hydration level could be insufficient for healthy living.

NORMAL HEALTHY RANGE OF BODY FAT & TOTAL BODY WATER PERCENTAGE

MALE			FEMALE		
AGE	BODY FAT	BODY WATER	AGE	BODY FAT	BODY WATER
18-29	8-19%	60%	18-29	20-30%	50%
30-39	11-21%		30-39	22-32%	
40-49	13-23%		40-49	24-34%	
50-59	15-25%	55%	50-59	26-36%	47%
60-69	17-27%		60-69	28-38%	
		50%			45%

Fig. 1

	BELOW NORMAL BODY WATER %	NORMAL & ABOVE BODY WATER %
LOW BODY FAT %	Aim for a normal BMI value, ensure sufficient water intake	Aim for a normal BMI value
HEALTHY BODY FAT %	Weight and body fat OK but insufficient water intake	Ideal balance body fat, muscle mass and weight
HIGH BODY FAT %	Change of diet required, ensure sufficient water intake	Change of diet required

Fig. 2

Body water measurement results are influenced by the proportion of body fat and muscle. If the proportion of body fat is high, or the proportion of muscle is low then the body water results will tend to be low.

It is important to remember that measurements such as body weight, body fat and body water are tools for you to use as part of your healthy lifestyle. As short term fluctuations are normal, we suggest you chart your progress over time, rather than focus on a single days reading.

This product is for domestic use only.

Professional medical guidance should always be sought before embarking on diet and exercise programmes.


This Salter Scale is calibrated for adult use (age 16+ years).

PREGNANT WOMEN SHOULD ONLY USE THE WEIGHT FUNCTION.

CAUTION: Do not use this product if you have a pacemaker or other medical device fitted. If in doubt, consult your doctor.

WHAT DOES MY BMI VALUE MEAN?

Body Mass Index (BMI) is an index of weight-for-height that is commonly used to classify underweight, overweight and obesity in adults. This scale calculates your BMI value for you. The BMI categories shown in the chart and table below are recognised by the WHO (World Health Organisation) and can be used to interpret your BMI value.

Your BMI category can be identified using the table below.

BMI Category	BMI Value Meaning	Health Risk Based Solely On BMI
< 18.5	Underweight	Moderate
18.6-24.9	Normal	Low
25-29.9	Overweight	Moderate
30+	Obese	High


IMPORTANT: If the chart/table shows that you are outside of a 'normal' healthy weight range, consult your doctor before taking any action.

BMI classifications are for adults (age 20+ years) only.

WEIGHT & BODY FAT READINGS

Position the scale on a flat level surface.

- 1) Press the **ON/OFF** button to turn the scale on.
- 2) Select your user number by pressing the **▲** or **▼** buttons.
- 3) **WAIT** while the display reconfirms your personal data then shows a zero reading.
- 4) When zero is displayed, **WITH BARE FEET**, step onto the platform and stand still. Ensure your feet are in good contact with the 2 metal pads on the platform.
- 5) After 2-3 seconds, your weight will be displayed followed by your body fat percentage, a total body water percentage, muscle mass percentage and BMI.

WARNING INDICATORS

'LO' = replace battery.

'Err2' = out of body fat range or poor feet contact.

'Err' = overload.


TECHNICAL SPECIFICATION

100 cm - 220 cm
3'3"¹/₂ - 7'2"¹/₂
16 - 80 years
Male/Female

10 user memory
d = 0.1 % body fat
d = 0.1 % body water
d = 0.1 % muscle mass


WEIGHT READING ONLY

Position the scale on a flat level surface.


- a) Tap the centre of the platform firmly to activate the scale.
- b) When '8' appears, remove foot from platform and **WAIT** for display to show a zero reading.
- c) When zero is displayed, step onto the platform and stand still.
- d) After 2-3 seconds, your weight will be displayed.

BATTERY


Requires 3 x AAA batteries.

WEIGHT MODE

Set the switch on underside of product to: KG, ST, or LB.

POSITION

Position scale on flat surface.

ADVICE FOR USE AND CARE

- Always weigh yourself on the same scale placed on the same floor surface. Do not compare weight readings from one scale to another as some differences will exist due to manufacturing tolerances.
- Placing your scale on a hard, even floor will ensure the greatest accuracy and consistency.
- It is recommended that you measure at the same time of the day, preferably early evening before a meal for the most consistent results.
- Your scale rounds up or down to the nearest increment. If you weigh yourself twice and get two different readings, your weight lies between the two.
- Clean your scale with a damp cloth. Do not use chemical cleaning agents.
- Do not allow your scale to become saturated with water as this can damage the electronics.
- Treat your scale with care - it is a precision instrument. Do not drop it or jump on it.

TROUBLESHOOTING GUIDE

- You must have bare feet when using body fat and total body water measurements.
- The condition of your skin on the bottom of your feet can affect the reading. To get the most accurate and consistent reading, wipe your feet with a damp cloth, leaving them slightly moist before stepping on the scale.
- When taking body fat readings ensure your feet are in good contact with the 2 metal pads on the platform.

GUARANTEE

This product is intended for domestic use only. Salter will repair or replace the product, or any part of this product, (excluding batteries) free of charge if within 15 years of the date of purchase, it can be shown to have failed through defective workmanship or materials. This guarantee covers working parts that affect the function of the scale. It does not cover cosmetic deterioration caused by fair wear and tear or damage caused by accident or misuse. Opening or taking apart the scale or its components will void the guarantee. Claims under guarantee must be supported by proof of purchase and be returned carriage paid to Salter (or local Salter appointed agent if outside the UK). Care should be taken in packing the scale so that it is not damaged while in transit. This undertaking is in addition to a consumer's statutory rights and does not affect those rights in any way. For UK Sales and Service contact HoMedics Group Ltd, PO Box 460, Tonbridge, Kent, TN9 9EW, UK. Helpline Tel No: (01732) 360783. Outside the UK contact your local Salter appointed agent.

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