

Advise for the best results by using the Ultra Slim

3 other key points regarding the use of the USL-1000 Ultra Slim

1. First is the pins and needle feeling inside whilst carrying out the treatment. This pins and needle feeling is the proof that the cavitation is working and that the fat is getting exploded. It's maybe not for the user to decide which kind of fat they have to set up the right parameters. So, THE PINS and NEEDLE FEELING should be the consequence that the right parameters have been chosen.
2. The second thing is the importance of taking water, before and after the treatment. The water is the fluid to dissolve the exploded fat, so the lymphatic system, has a better chance to get the fat out of the body. Sometimes, the user takes lots of water just before the treatment and, this is not good if the treatment is going to be made on the abdomen area. The water should be taken, at least, 45 to 1 hour before the treatment and, during the 24 to 48 hours after the treatment.
3. One more important thing is that the user, after having a treatment and if he chose the right parameters, he will see that his urine, a couple of hours later, will be a bit darker than usual. This is not a problem, it is another proof than the cavitation treatment, really exploded the fat and the lymphatic system is taking it out of the body.

The most important thing though to getting the best results from the product, is drinking plenty of water at least 1 hour before the treatment and then after the treatment for the next couple of days, also to carry out some light exercise after the treatment. This will all help to improve the results, but the thing to remember is this is not a miracle cure and the time taken to see results will vary from 1 person to the next. The key is to persevere to give the product more time and to ensure that they are drinking plenty of water, eating healthily and doing a moderate amount of exercise, which will speed up the results.